

Mobilscooter priser 2023				
Product	Part number	EAN Code	Product Photo	Vejl. Pris Inkl moms

ATTO
Mobility Scooter - Classic White
upto 120 KG User Weight / Speed upto 6.4 KM
Including Standard Flight Battery 249.6 WH

Contains 1 product	
ATTO	
AT01-101-B2-0	7290110012878



28.995,-

ATTO SPORT
Mobility Scooter - Dark Grey
upto 120 KG User Weight / Speed upto 10 KM
Including Standard Flight Battery 249.6 WH

Contains 1 product	
ATTO SPORT	
AT02-100-B2-0	7290110013110



34.995,-

ATTO SPORT MAX
Mobility Scooter - Dark Grey
upto 136 KG User Weight / Speed upto 6.4 KM
Including Standard Flight Battery 249.6 WH

Contains 2 products	
ATTO SPORT	
AT02-100-B2-0	7290110013110
MAX KIT	
600-004788	7290110012854



37.995,-

NOTE: MAXI KIT NEEDS TO BE INSTALLED BY PARTNER



ATTO Family Comparison Chart


ATTO • ATTO SPORT • ATTO SPORT MAX • SHABB ATTO • SHABB ATTO MAX



FEATURE	ATTO			ATTO SPORT			ATTO SPORT MAX			SHABBATTO			SHABBATTO MAX		
Max Speed	6.4 km/h / 4 mph			10 km/h / 6.2 mph			6.4 km/h / 4 mph			Regular use: 10 km/h / 6.2 mph During Shabbat: 6.4 km/h 4 mph			Regular use: 6.4 km/h / 4 mph During Shabbat: 6.4 km/h 4 mph		
Max User Weight	120 kg / 265 lb			120 kg / 265 lb			136 kg / 300 lb			120 kg / 265 lb			136 kg / 300 lb		
Battery Size	249.6 wh			249.6 wh			249.6 wh			249.6 wh			249.6 wh		
Driving Distance Per Charge	Up to 20 km / 12.5 miles			Up to 20 km / 12.5 miles			Up to 17 km / 10.5 miles			Regular use: up to 20 km / 12.5 miles During Shabbat: up to 13 km / 8 miles			Regular use: up to 17 km / 10.5 miles During Shabbat: up to 10.5 km / 6.5 miles		
Dimensions	Drive mode: L 47" / H 35" / W 22" Trolley mode: L 15.3" / H 28" / W 16.5"			Drive mode: L 47" / H 35" / W 22" Trolley mode: L 15.3" / H 28" / W 16.5"			Drive mode: L 47" / H 35" / W 22" Trolley mode: L 15.3" / H 28" / W 16.5"			Drive mode: L 47" / H 35" / W 22" Trolley mode: L 15.3" / H 28" / W 16.5"			Drive mode: L 47" / H 35" / W 22" Trolley mode: L 15.3" / H 28" / W 16.5"		
Cruise Control (Throttle Not Needed)	No			Yes			Yes			Yes			Yes		
Automatic Slowdown	Slopes			Slopes / Turns			Slopes			Slopes / Turns			Slopes		
Front Disk Brake for Extra Safety	No			Yes			Yes			Yes			Yes		
Display	LCD / Programmable			LCD / Programmable			LCD / Programmable			LCD / Programmable			LCD / Programmable		
Tires	Airless PU, flat free tires			Airless PU, shock absorbing, flat free tires			Airless PU, heavy duty, dense, shock-absorbing, flat-free tires			Airless PU, shock absorbing, flat free tires			Airless PU, heavy duty, dense, shock-absorbing, flat-free tires		
Antitheft Wireless Locking FOB Key	No			Yes			Yes			Yes			Yes		
Embedded Front Headlights	No			Yes			Yes			Yes			Yes		
Seat and Backrest / Grips	PU Seat / Grips			PU Leather Seat / Grips			PU Leather Seat / Grips			PU Leather Seat / Grips			PU Leather Seat / Grips		
Rear Light Lifting Handle	No			Yes			Yes			Yes			Yes		
Eco Mode	No			Limit speed / Increase range			No			Limit speed / Increase range			No		
Left-hand Option	Optional			Included			Included			Included			Included		
Shabbat Remote	No			No			No			Yes			Yes		
Weight	Front Half: 26.5 lb / 12 kg Rear Half: 37.5 lb / 17 kg Assembled: 64 lb / 29 kg Battery: 4 lb / 1.8 kg			Front Half: 27.5 lb / 12.5 kg Rear Half: 38.5 lb / 17.5 kg Assembled: 66 lb / 30 kg Battery: 4 lb / 1.8 kg			Front Half: 27.5 lb / 12.5 kg Rear Half: 38.5 lb / 17.5 kg Assembled: 66 lb / 30 kg Battery: 4 lb / 1.8 kg			Front Half: 27.5 lb / 12.5 kg Rear Half: 38.5 lb / 17.5 kg Assembled: 66 lb / 30 kg Battery: 4 lb / 1.8 kg			Front Half: 27.5 lb / 12.5 kg Rear Half: 38.5 lb / 17.5 kg Assembled: 66 lb / 30 kg Battery: 4 lb / 1.8 kg		
Dimensions trolley Mode	H 28"	L 15.3"	W 16.5"	H 30.1"	L 15.3"	W 17.5"	H 30.1"	L 15.3"	W 17.5"	H 30.1"	L 15.3"	W 17.5"	H 30.1"	L 15.3"	W 17.5"
	H 72 cm	L 39 cm	W 42 cm	H 77.5 cm	L 39 cm	W 44.5 cm	H 77.5 cm	L 39 cm	W 44.5 cm	H 77.5 cm	L 39 cm	W 44.5 cm	H 77.5 cm	L 39 cm	W 44.5 cm
Dimensions Drive Mode	H 35"	L 15.3"	W 22"	H 35"	L 15.3"	W 22.8"	H 35"	L 15.3"	W 22.8"	H 35"	L 15.3"	W 22.8"	H 35"	L 15.3"	W 22.8"
	H 90 cm	L 39 cm	W 56 cm	H 90 cm	L 39 cm	W 58 cm	H 90 cm	L 39 cm	W 58 cm	H 90 cm	L 39 cm	W 58 cm	H 90 cm	L 39 cm	W 58 cm

ATTO Armrests (Classic White)		600-004225	7290110010034	■ ATTO □ ATTO_{SPORT} □ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	3825,-
ATTO SPORT Armrests (Dark Grey)		600-104225	7290110012038	□ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	3825,-
Seat Cushion - Black		600-004230	7290110011413	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	995,-
Seat Cushion - Blue		600-004232	7290110011482	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	995,-
Seat Cushion - Red		600-004231	7290110011475	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	995,-
Seat Cushion Leather		600-104837	7290110013318	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	1495,-
Luxury Orthopedic Backrest Leather		600-004841	7290110013356	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	€ 2.495
Pet Carrier & Cushion Black		600-004840	7290110013349	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	1495,-
Luxury Swivel Seat Leather		600-004842	7290110013363	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	1495,-
Seat & Support Belt		600-004839	7290110013332	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	800,-
Oxygen Bottle Holder & Case		600-004843	7290110013370	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	995,-
Essentials Pouch		600-004235	7290110011444	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	1095,-
Front / Rear Basket Bag		600-004236	7290110011499	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	1695,-
Carryal & Cushion		600-004233	7290110011451	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	1495,-
Carryall & Cushion Leather		600-104838	7290110013325	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	1995,-
Flight Kit		600-004234	7290110011468	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	2495,-
Cane / Crutches Holders		600-004222	7290110010041	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	649,-
Mobile Phone Holder		600-004229	7290110011505	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	€ 475
ATTO Left Hand Throttle Kit		600-004395	7290110011611	■ ATTO □ ATTO_{SPORT} □ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	1795,-
Standard Flight Battery - 249.6 Wh *Standard on ATTO / ATTO SPORT		110-001192-2P-A	7290110010010	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	5250,-
Cable Kit		600-004238	7290110011437	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	495,-
Accessories Knobs Long		600-004733	7290110012250	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	195,-
Rear Light Lifting Handle *Standard on ATTO SPORT / SHABBATTO + MAX Series		600-104609	7290110012267	■ ATTO ■ ATTO_{SPORT} ■ ATTO_{MAX} ■ ATTO_{SPORT} ■ ATTO_{MAX}	495,-


NOTE: For use of multiple accessories



SPARE PARTS - ATTO

Item	Photo	Description	Part number	EAN Code	Contains	Vejl. Pris Inkl moms
1		ATTO / ATTO SPORT Telescope Wire	700-003993	7290110011147	Telescope Wire X1 Small Zip-Tie X1 Medium Zip-Tie X2	400,-
2		ATTO Telescope Assembly Blue	700-003996	7290110011178	Telescope Assembly X1	1450,-
3		ATTO Telescope Assembly Gray	700-003996-G	7290110012229	Telescope Assembly X1	1450,-
4		ATTO Back Shell	700-004241	7290110011253	HB Rear Shell Over Mold Assembly X1 HB USB Cover X1PT Cross Pin Head Screw X12	500,-
5		ATTO HB-Grip R-L	700-004281	7290110011322	Spring Connector X2 HB Grip Machining X2 HB Grip L X1 HB Grip R X1 Spring Pin - A2 X2	1050,-
6		ATTO Front Shell	700-004279	7290110011314	HB Front Shell Over Mold Assembly X1 HB Front Logo X1 PT S1 20x2 Cross Pin Head X2 PT Cross Pin Head Screw X5	400,-
7		ATTO / ATTO SPORT Handle Axis Lock Set	700-004269	7290110011307	HB Handle Axis Pin X1 Handle Axis Pin Lock X1 Hex. Socket Set Screw X1	200,-
8		ATTO / ATTO SPORT Separation ASM (one unit)	700-004240	7290110011246	Separation Assembly X1	1500,-
9		ATTO Folding Handle Mechanism	700-004282	7290110011338	Spring Connector X2 HB Chassis Tube Retainer X2 Spring Pin - A2 X2 Handles Retaining Spring X1	150,-
10		ATTO Front Right Bay	700-004286	7290110011376	Right Bay X1 Bumper Right X1 Right Cover X1 USB Lid X1 Cross Pin Head Screw X4 PT Cross Pin Head Screw X5 Black Coat Cross Pin X2 Right Bay Label-English X2	600,-
11		ATTO / ATTO SPORT Chassis Lock Plunger	700-004289	7290110011406	Plunger Casting X1 Chassis Lock Plunger X1 Lee Spring UC202B0055 X1	150,-
12		ATTO Handle Bar	700-004181	7290110011208	Handle Bar Assembly For Assembly Line X1 Handle Bar Label-English X1	5000,-
13		ATTO Handle Bar (With LCD Screen)	700-004815	7290110013097	Handle Bar Assembly For Assembly Line X1 Handle Bar Label-English X1	5000,-
14		ATTO Battery Bay	700-004028	7290110011185	Front Left Assembly X1 Battery Bay Label-English X1 Trolley Label L-English X1	700,-
15		ATTO Throttle	700-003992	7290110011130	Throttle Harness Right X1	400,-
16		ATTO Controller Rev05 (For LCD Screen)	700-004871	7290110013455	Controller Box ASM - Rev C X1 Cross Countersunk Head Screw X3 Zip-Tie - Small X5	5000,-
17		ATTO Front Wheel Motor (New Version)	700-004819	7290110013257	Limx Hub Motor Rev A X1 Tab Washer for M12 X1 Hex Nut - A2 Black Coating X1 Torque Nut - A2 Black Coating X1 Hex Socket Cap Head - Black X2 Zip-Tie - Small X4 Motor Label X1	7800,-
18		ATTO / ATTO SPORT T-Plunger Kit	700-003989	7290110010096	T-Plunger X1	300,-
19		ATTO Backrest	700-003988	7290110010089	Backrest X1 Half Tooth With Glue Screw PM 3.0*8mmxglue 4mm X4 Half Tooth With Glue Screw PM 3.0*30mmxglue 4mm X4	1500,-
20		ATTO Armrest Dummy R-L	700-003991	7290110011123	Armrest Dummy L X1 Armrest Dummy R X1	500,-
21		ATTO Chair Assy Kit	700-003995	7290110011161	Chair Assembly X1	9500,-
22		ATTO Wheel and Plastic Cover	700-004245	7290110011291	Wheel And Bearing Assy X1 Wheel Axis Plastic Cover X1	1000,-

21		ATTO Wheel Arm ASM	700-004244	7290110011284	Wheel Housing Extrusion Weld Assy X1 Wheel Bearing Cap X1 Wheel Securing Nut X1 Torque Nut - A2 X1 Hex socket Countersunk head X1	900,-
22		ATTO Chair Dress	700-004243	7290110011277	Dress X1 PT Cross Pan Head Screw X6	300,-
23		ATTO Chair Back Cover	700-004242	7290110011260	Backcover Top Assembly X1 PT Cross Pan Head Screw X6 PT ST Zinc Cross Pan Head X2 Flat Washer D7x6.3 27x10.5mm X8	1500,-
24		ATTO Left Slider Blue	700-004231	7290110011239	Side Face For Assy X1 Slider Plastic Terminal X1 Hinge Roller Assy X1 Lee Spring LCG28B0055 X2 Hex Socket Countersunk Head X2 Hex Socket Cap Head - Black X1	950,-
		ATTO Left Slider Gray	700-004231-G	7290110012199	Side Face For Assy X1 Slider Plastic Terminal X1 Hinge Roller Assy X1 Lee Spring LCG28B0055 X2 Hex Socket Countersunk Head X2 Hex Socket Cap Head - Black X1	950,-
25		ATTO Rear Profile Hole Cap	700-004283	7290110011345	Rear Profile Hole Cap X5	120,-
26		ATTO Foot Support - Bottom Cover	700-004284	7290110011352	Foot Support X1 Left Bottom Cover X1 Right Bottom Cover X1 PT Cross Pan Head Screw X10 Flat Washer X2	1000,-
27		ATTO L-Profile Top Cap	700-004285	7290110011369	L Profile Top Cap X1	100,-
28		ATTO / ATTO SPORT Hinge Roller	700-004287	7290110011383	Hinge Roller Assy X1	150,-
29		ATTO Armrest Lock R-L	700-004288	7290110011390	Armrest Release Button L Assy X1 Armrest Release Button R Assy X1 Black Coat Cross Pan h 3x4	250,-
30		ATTO Right Slider Blue	700-004224	7290110011222	Right Side Face Assy X1 Hinge Roller Assy X1 Hex Socket Cap Head - Black X1 Hex Socket Countersunk Head 5 X2 Slider Plastic terminal	600,-
		ATTO Right Slider Gray	700-004224-G	7290110012205	Right Side Face Assy X1 Hinge Roller Assy X1 Hex Socket Cap Head - Black X1 Hex Socket Countersunk Head 5 X2 Slider Plastic terminal	600,-
31		ATTO Chassis Releaser	700-004029	7290110011192	Chassis Releaser Assy X1 Front Chassis Locking Pin Stopper X2 Front Chassis Locking Pin X1 Lee Spring LCG30D155M X1 Spring Pin - A2 X2	450,-
32		ATTO / ATTO SPORT Rear Handle (Tail)	700-003987	7290110010072	Tail Handle X1 PT Cross Pan Head Screw X2 PT ST Black Coat Cross Pan h 3x4	300,-
33		ATTO Seat Upholstery	700-003986	7290110010065	Seatrest X1 PT A2 Black Coat X6	1300,-
34		ATTO / ATTO SPORT Chair Frame	700-004298	7290110011420	Steel Chair Frame X1 Seat Support Hinge X2 Flange Bearing IGUS GPM-5517-04 X4 Hex Socket Pan Head Screw X2 A Flat Washer X2	1700,-
35		ATTO Rod Latch Button	700-004188	7290110011215	Rod Latch Button Frame X1 Rod Latch Button Cover X1 PT ST Black Coat Cross Pan Head X1 A2 Cross Pan Head Screw X2 Lee Spring LCG24C295 X1	200,-
36		ATTO Rear Assy	700-004299	7290110011604	Rear Assembly X1	19900,-
37		ATTO Front Assy (without Battery)	700-004300	7290110011659	Front Assembly Without Battery	23500,-
38		ATTO Top Buttons Module	700-003997	7290110011673	HB Top Buttons Module Assembly X1 PT ST Zinc Cross Pan Head X2	400,-
No Number		ATTO / ATTO SPORT Universal Charger	700-003835	7290110010003	Universal Charger Not included Plug + Cable	995,-
40		ATTO Rear Lifting Handle	700-004226	7290110010027	Lifting Handle X1 Hanging knob X2 Instructions Lifting Handle Accessory X1	300,-
41		ATTO HB Grip Right+Left	700-003434	7290110012106	HB Grip R HB Grip L	150,-
42		ATTO Chair Joint	700-003120	7290110012113	Chair Joint (Assembled w/ 2 rod bearings)	800,-
43		ATTO / ATTO SPORT Slider Plastic Terminal	700-003316	7290110011826	Slider Plastic terminal	200,-
44		ATTO Display module	700-003784	7290110012120	HB Display Module Assembly	1400,-
45		ATTO / ATTO SPORT Accessories Knobs Short	700-004686	7290110012137	Accessories knob X2	50,-
46		ATTO Telescope Locking Lever Blue	700-004419	7290110011864	Telescope Handle Holding Pin X1 Telescope Handle Casing X1 Telescope Latch Plastic Teeth X1 Telescope Latch Handle Rubber X1 Telescope Latch Tongue Spring X1	300,-
		ATTO Telescope Locking Mechanism Gray	700-004419-G	7290110012212	Telescope Handle Holding Pin X1 Telescope Handle Casing X1 Telescope Latch Plastic Teeth X1 Telescope Latch Handle Rubber X1 Telescope Latch Tongue Spring X1	300,-
47		ATTO / ATTO SPORT Chair Handle Kit	700-004420	7290110011871	Chair Handle Assembly	350,-
48		ATTO Steering Plastic Covers Kit	700-004000	7290110011765	Apple Cover Front Apple Cover Back	150,-
50		ATTO / ATTO SPORT Armrest Installation Pin	700-004325	7290110011840	Armrest Installation Pin	300,-
51		ATTO Rod Link Set	700-004782	7290110012822	Link R for Assembly Line X1 Link L for Assembly Line X0	400,-



SPARE PARTS - ATTO SPORT

Item	Photo	Description	Part number	EAN Code	Content	Vejl. Pris Inkl moms
1		ATTO / ATTO SPORT Telescope Wire	700-003993	7290110011147	Telescope Wire X1 Small Zip Tie X3 Medium Zip Tie X2	400,-
2		ATTO SPORT Telescope Assembly	700-103996	7290110012274	Telescope Assembly X1	1450,-
3		ATTO SPORT Back Shell	700-104241	7290110012281	HB Rear Shell Over Mold Assembly X1 HB USB Cover X1 PT Cross Pin Head Screw X12	600,-
4		ATTO SPORT Front Shell	700-104279	7290110012298	HB Front Shell Over Mold Assembly X1 HB Front Logo X1 PT ST 20x20 Cross Pin Head X2 PT Cross Pin Head Screw X8	425,-
5		ATTO / ATTO SPORT Handle Axis Lock Set	700-004269	7290110011307	HB Handle Axis Pin X1 Handle Axis Pin Lock X1 A2 Hex Socket Set Screw X1	200,-
6		ATTO / ATTO SPORT Separation ASM (one unit)	700-004240	7290110011246	Separation Assembly X1	1500,-
7		ATTO SPORT Folding Handle Mechanism	700-104282	7290110012304	Spring Connector X2 HB Chassis Tube Retainer X2 Spring Pin - A2 X2 Handles Retaining Spring X1	150,-
8		ATTO SPORT Front Right Bay	700-104286	7290110012311	Right Bay X1 Front Led Assy R Right Cover X1 USB Lid X1 A2 Cross Pin Head Screw X4 PT Cross Pin Head Screw X5 A2 BLACK Coat Cross Pin X2	675,-
9		ATTO / ATTO SPORT Chassis Lock Plunger	700-004289	7290110011406	Plunger Casing X1 Chassis Lock Plunger X1 Lee Spring LC020B0055 X1	150,-
10		ATTO SPORT Handle Bar (L+R)	700-104834	7290110012328	Handle Bar Assembly For Assembly Line X1 Handle Bar Label English X1	6000,-
11		ATTO SPORT Battery Bay	700-104028	7290110012335	Front Left Assembly X1 Battery Bay Label English X1 Trolley Label L English X1	950,-
No Number		ATTO / ATTO SPORT Universal Charger	700-003835	7290110010003	Universal Charger Batter Not included Plug + Cable	995,-
13		ATTO SPORT Throttle	700-103992	7290110012359	ATTO SPORT Throttle Assembly X1	350,-
14		ATTO SPORT Controller Rev05	700-104870	7290110013448	ATTO SPORT Controller Assembly X1 A2 Cross Pin Head screw X2	4000,-
15		ATTO SPORT Armrest Dummy R-L	700-103991	7290110012373	Armrest Dummy L X1 Armrest Dummy R X1	500,-
16		ATTO SPORT Chair Assy Kit	700-103995	7290110012380	ATTO SPORT Chair Assembly X1	16500,-
17		ATTO Sport Wheels assembly Left+Right	700-104245	7290110012397	ATTO SPORT Wheel R X3 ATTO SPORT Wheel L X3	1700,-
18		ATTO SPORT Wheel Arm ASM	700-104244	7290110012403	Wheel Housing Extrusion Weld Assy X1 Wheel Bearing Cap X1 Wheel Securing Nut X1 Torque Nut - A2 X1 A2 Hex Socket Countersunk Head X1	900,-
19		ATTO SPORT Chair Dress	700-104243	7290110012410	ATTO SPORT Dress X1 PT Cross Pin Head Screw X6	300,-
20		ATTO SPORT Chair Back Cover	700-104242	7290110012427	Backcover Top Assembly X1 PT Cross Pin Head Screw X6 PT ST 20x20 Cross Pin Head X2	1450,-
21		ATTO SPORT Left Slider	700-104231	7290110012434	Side Face For Assy X1 Slider Plastic Terminal X1 Hinge Roller Assy X1 Lee Spring LC020B0055 X1 A2 Hex Socket Countersunk Head X2 A2 Hex Socket Cap Head - Black X1	950,-
22		ATTO SPORT Rear Profile Hole Cap	700-104283	7290110012441	Rear Profile Hole Cap X5	120,-

23		ATTO SPORT Foot Support - Bottom Cover	700-104284	7290110012458	Foot Support X1 Left Bottom Cover X1 Right Bottom Cover X1 PT Cross Pin Head Screw X30 A2 Flat Washer X2	1000,-
24		ATTO SPORT L-Profile Top Cap	700-104285	7290110012465	L Profile Top Cap X1	100,-
25		ATTO / ATTO SPORT Hinge Roller	700-004287	7290110011383	Hinge Roller Assy X1	150,-
26		ATTO SPORT Armrest Lock R-L	700-104288	7290110012472	Armrest Release Button L Assy X1 Armrest Release Button R Assy X1 PT ST. Black Coat Cross Pan H X4	250,-
27		ATTO SPORT Right Slider	700-104224	7290110012489	Right Side Face Assy X1 Slider Plastic Terminal X1 Hinge Roller Assy X1 A2 Hex Socket Cap Head- Black c. X1 A2 Hex Socket Countersunk Head S X2	600,-
28		ATTO SPORT Chassis Releaser	700-104029	7290110012496	Chassis Releaser Assy X1 Front Chassis Locking Pin Stopper X2 Front Chassis Locking Pin X1 Lee Spring LCG3803SM X1 Spring Pin -A2 X2	450,-
29		ATTO / ATTO SPORT Rear Handle (Tail)	700-003987	7290110010072	Tail Handle X1 PT Cross Pin Head Screw X2 PT ST. Black Coat Cross Pan H X4	300,-
30		ATTO / ATTO SPORT T-Punger Kit	700-003989	7290110010096	T-Punger X1	300,-
31		ATTO SPORT Front Wheel Motor	700-104836	7290110013301	ATTO SPORT Front Wheel Motor X1 Hub Motor NPT X1 Motor Spacer X1 Tab Washer For M12 X3 Hex Nut - A2 Black Coating X1 Helical Spring Lock Washer B 12 X1	7300,-
32		ATTO SPORT Backrest	700-103988	7290110012519	Backrest X1 Half Tooth With Glue Screw PM 3.0*8mm 4mm X4 Half Tooth With Glue Screw PM 3.0*10mm 4mm X4	1500,-
33		ATTO SPORT Seat Upholstery	700-103986	7290110012526	Seatrest X1 PT A2 Black Coat HS	1300,-
34		ATTO / ATTO SPORT Chair frame	700-004298	7290110011420	Steel Chair Frame X1 Seat Support Hinge X2 Flange Bearing G615 GFM J517.04 X4 A2 Hex Socket Pan Head Screw X2 Flat Washer X2	1700,-
35		ATTO SPORT Rear Assy	700-104299	7290110012533	Rear Assembly X1	19900,-
36		ATTO SPORT Assy (Without Battery)	700-104300	7290110012540	Front Assembly Without Battery X1	23500,-
37		ATTO SPORT Top Buttons Module	700-103997	7290110012557	HB Top Buttons Module Assembly X1 PT ST. Zinc Cross Pan Head X2	400,-
38		ATTO SPORT Handle grip	700-104748	7290110012564	Plastic Grip With Leather Cover X1 Grip Handle Cup X1 Grip Handle Locker X1	250,-
39		ATTO SPORT Handle brake unit	700-104749	7290110012571	Brake Handle X1 Lever Screw X1 Brake Handle Locking Nut X1 Brake Cable Spring X1 Brake Spring Protector X1 Brake Cable + Cover X1 Brake Bundle Sleeve X1 Brake Cable Shrink X2	400,-
40		ATTO / ATTO SPORT Armrest Installation Pin	700-004325	7290110011840	Armrest Installation Pin	300,-
41		ATTO SPORT Display module	700-103784-R	7290110012588	HB Display Module Assembly X1	1300,-
42		ATTO SPORT Chair joint	700-103120	7290110012595	Chair Joint (Assembled with 2 rod bearings) X1	800,-
43		ATTO / ATTO SPORT Slider Plastic Terminal	700-003316	7290110011826	Slider Plastic Terminal	200,-
44		ATTO SPORT Direction Button Assy	700-104750	7290110012601	Direction Button Assy X1	450,-
45		ATTO SPORT Caliper	700-104751	7290110012618	Caliper X1	300,-
46		ATTO SPORT Disk Brake	700-104756	7290110012625	ATTO SPORT Front Disk X1 Screw 92137M414 M5 X30 Black Oxide X5	200,-
47		ATTO SPORT Telescope Locking Lever	700-104419	7290110012632	Telescope Handle Holding Pin Telescope Handle casing Telescope Latch Plastic Teeth Telescope Latch Handle Rubber Telescope Latch Torque Spring	300,-
48		ATTO SPORT Steering Plastic Covers Kit	700-104000	7290110012847	Apple Cover Front Apple Cover Back	250,-
49		ATTO / ATTO SPORT Chair handle kit	700-004420	7290110011871	Chair Handle Assembly	350,-
50		ATTO / ATTO SPORT Accessories Knobs Short	700-004686	7290110012137	Accessories Knob X2	50,-
51		ATTO SPORT Grip Spacer Rubber Kit	700-104753	7290110012649	ATTO SPORT Grip Spacer Rubber Left X1 ATTO SPORT Grip Spacer Rubber Right X1	100,-
52		ATTO SPORT Handle Bar Rubber Kit	700-104754	7290110012656	Brake Cable Holder X1 Brake Spring Protector X1 Brake Handle Spacer X1	100,-
53		ATTO SPORT Rod Link Set	700-104782	7290110012854	Link R For Assembly Line X1 Link L For Assembly Line X1	400,-
54		ATTO SPORT Front Led kit	700-104755	7290110012663	Front Led Assy R X1 Front Led Assy L X1	400,-
55		ATTO SPORT Rod Latch Button	700-103994	7290110012342	Rod Latch Button Frame X1 Rod Latch Button Cover X1 PT ST. Black Coat Cross Pan Head X1 A2 Cross Pin Head Screw X2 Lee Spring LCG3AC095 X1	200,-

The purchase of Units, Accessories, Spare Parts and Promotion Materials are Ex. Works Warehouse Europe / Warehouse UK (United Kingdom) and does not include Shipping, VAT, Customs, Duty, Handling and / or labour costs.

The range within the scooter operates on a single battery charge can vary greatly depending on various factors, such as: Load Capacity (User weight + Carry Weight), Weather, Going Uphill and frequent speed changes.

The prices of end-users may differ per country, as each country has a different VAT rate.

*Movinglife has a strict policy regarding our recommended prices for end-users.
Price dumping on the internet is not appreciated by Movinglife.*

All right reserved, MovingLife Ltd / MovingLife UK Ltd. © 2023



ATTO Drive
Will take you anywhere
you want to go



ATTO Trolley
Folds easily and
wheels like a trolley



ATTO Compact
Opens to full size
in 10 seconds



ATTO Split
Splits into 2 parts,
for easy lifting and storing

movinglife®
life changing solutions